



N I R L A B

ON-DEMAND SOLUTIONS

The Revolution of Substance Analysis

Empowering your organization to develop mobile lab applications that enable instant analysis of anything, anywhere.

THE COMPANY

NIRLAB AG

Mobile AI lab to analyze anything anywhere instantly.

Established in 2018, NIRLAB AG, a Swiss spin-off from the University of Lausanne, has revolutionized the way professionals and organizations analyze materials using NIR and Raman spectroscopy and advanced machine learning.


"With our digital ecosystem we bring high precision labs to the field and enable rapid decision making based on trustable data."

Florentin Coppey, Founder NIRLAB AG



GLOBAL PRESENCE

Scanning Substances in +35 Countries across the Globe



Pioneering Scientific Innovation

Our esteemed partnership with the **Forensic Institute of the University of Lausanne** in Switzerland has made us a recognized name in global scientific circles.

Our **contributions to top-tier forensic, science, and pharmaceutical journals** validate our commitment to advancing knowledge and pushing technological frontiers.

Unil

UNIL | Université de Lausanne

Florentin Coppey,^a Cédric Schelling,^{b, c} Jean-Luc Veuthey,^{b, c} and Pierre Esseiva^{a*}

^a University of Lausanne, School of Criminal Justice, CH-1015, Lausanne, Switzerland, e-mail: pierre.esseiva@unil.ch

^a University of Lausanne, School of Criminal Justice, CH-1015, Lausanne, Switzerland,
e-mail: pierre.esseiva@unil.ch

Robert Deschenaux on the occasion of his retirement

published by Wiley-VHCA AG. This is an open access article under the terms of the [Creative Commons Attribution License](#), which permits use and distribution in any medium, provided the original work is properly cited. NoDerivs License, which permits use and distribution in any medium, provided the original work is properly cited and no modifications or adaptations are made.

been significant interest within the forensic community regarding the use of portable NIR devices to provide real-time results. This article introduces an innovative technology that integrates a handheld device, specifically, *Viavi MicroNIR*, with a server that encompasses a server responsible for data processing and a mobile

Journal of Pharmaceutical and Biomedical Analysis 202 (2021) 11415

Contents lists available at ScienceDirect

Journal of Pharmaceutical and Biomedical Analysis

journal homepage: www.elsevier.com/locate/jpba

New perspective for the in-field analysis of cannabis samples using handheld near-infrared spectroscopy: A case study focusing on the determination of Δ^9 -tetrahydrocannabinol

Riccardo Deidda^{a,*,1}, Florentin Coppey^{b,1}, Dhouha Damerji^{c,d}, Cédric Schelling^{c,d},
Lauren Coïc^a, Jean-Luc Veuthey^{c,d}, Pierre-Yves Sacré^a, Charlotte De Bleye^a,
Philippe Hubert^a, Pierre Esseiva^b, Éric Ziemons^a

² University of Liège (Université de Liège), CRM, Vibra-Santé HUB, Laboratory of Pharmaceutical Analytical Chemistry, B36 Tower 4 Avenue Hippocrate 15, 4000, Liège, Belgium

ARTICLE INFO

Article history:
Received 15 March 2021
Received in revised form 6 May 2021
Accepted 10 May 2021
Available online 19 May 2021

Keywords:
Cannabinoids
THC
Cannabis
Handheld
Near infrared spectroscopy
NIR

ABSTRACT

The aim of the present study was to explore the feasibility of applying near-infrared (NIR) spectroscopy to the quantitative analysis of Δ^9 -tetrahydrocannabinol (THC) in cannabis products using handheld devices. A preliminary study was conducted on different physical forms (entire, ground and sieved) on cannabis inflorescences in order to evaluate the impact of sample homogeneity on THC content predictions. Since entire cannabis inflorescences represent the most common types of samples found in both the pharmaceutical and illicit markets, they have been considered priority analytical targets. Two handheld NIR spectrometers were used to obtain NIR spectra of the samples. The effect of sample size on their predictive performance was compared. Six partial least square (PLS) models based on reference data obtained by HPLC-UV were built. The importance of the technical features of the spectrophotometers for quantitative applications was highlighted. The mid-cost system outperformed the low-cost system in terms of predictive performance, especially when analyzing entire cannabis inflorescences. In contrast,

The mid-cost system was selected as the best-suited spectrophotometer for this application. The number of cannabis inflorescence samples was augmented with new real samples, and a chemometric model based on machine learning ensemble algorithms was developed to predict the concentration of THC in 75 samples. Good predictive performance was obtained with a root mean squared error of prediction of 1.75% (w/w). The proposed method was used to compare the NIR predictions to the quantitative results obtained by HPLC-UV and evaluate the degree of accordance between the two analytical techniques. Each result fell within the established limits of agreement, demonstrating the feasibility of this chemometric model for analytical purposes.

Finally, resin samples were investigated by both NIR devices. Two PLS models were built by using a sample set of 45 samples. When the analytical performances were compared, the mid-cost spectrophotometer significantly outperformed the low-cost device for prediction accuracy and reproducibility.

1. Introduction

The analysis of cannabis samples mainly concerns two general areas: quality control laboratories (often for medicinal cannabis) and forensic laboratories (seized cannabis samples). The simplest medicinal cannabis samples available on the market consist of dried flower tips with the aim for use in various therapeutic indications (from multiple sclerosis to epilepsy) [1]. This product is

* Corresponding author.
E-mail address: riccardo.deidda@uliege.be (R. Deidda).

<https://doi.org/10.1016/j.jpba.2021.114150>
0731-7085/© 2021 Elsevier B.V. All rights reserved.

<https://doi.org/10.1016/j.jpba.2021.114150>
0731-7085/© 2021 Elsevier B.V. All rights reserved.

<https://doi.org/10.1016/j.jpba.2021.114150>
0731-7085/© 2021 Elsevier B.V. All rights reserved.

<https://doi.org/10.1016/j.jpba.2021.114150>
0731-7085/© 2021 Elsevier B.V. All rights reserved.

The background features a dark, almost black, space. A bright red, glowing ring or torus is positioned in the upper left, tilted at an angle. A beam of red light emanates from the center of this ring, extending downwards and slightly to the right. The light has a soft, ethereal quality with some internal texture. In the upper right, there are faint, scattered red dots, suggesting distant stars or particles.

The Solution for Your Organization

FIELDLAB for Onsite Analysis Tailored to your Needs

We provide your organization with high precision labs for the field in the fastest possible way.




Substance Library

The substance library consists of highly precise data models for the specific application areas and acts as reference library.




NIRLight

Substances can be easily scanned with NIR handheld devices. Device can be rapidly connected to a mobile NIRApp via Bluetooth.



NIRApp

NIRApp (iOS/Android) displays the analysis results instantly and enables standardized data collection and processing in the field.



NIRWeb

Desktop app and browser-accessible platform for data, license and application management. Various dashboards enable valuable insights across all scans and devices.




NIRCloud

Secure, high-speed server where analysis results and prediction models are trained and stored. Hosted on the University of Lausanne campus secured data center.

The Ecosystem for a Fast Go-to-Market

We provide an end-to-end solution to industries where instant identification and quantification of materials is required.



1 Data collection



R&D support
Data collection
Creating reference library
Associate laboratory results



2 Data Modelling & Calibration



Chemometric Suite
Modelling & dataset development
Data management & machine learning
Validation




3 Release & Go-To-Market

Fieldlab
Handheld NIR-device
Mobile app with instant reporting
Desktop app- and browser-accessible platform

How it works

From data acquisition to production, a highly automatized and customizable solution.





FIELDLAB by NIRLAB

High-precision handheld screening device for identifying and quantifying of substances across various industries.

NIRLight

THE HARDWARE


- + **Signal to Noise Ratio**
Among the highest in the field of handheld devices.
- + **Wireless, compact, rugged and ergonomic.**
Designed for use in the field as well as in the laboratory.
- + **IP65 and IP67 rated**
Made for wet and dusty environments.
- + **Destruction-free analysis**
Little or no sample preparation is needed. No special training required.
- + **No maintenance**
No maintenance is required. The glass and lamp are replaceable if broken.
- + **Long battery life**
USB charging and 10 hours of continuous use.
- + **Bluetooth and USB**
Simple and fast connectivity to tablet or PC.



NIRLab Mobile App

THE SOFTWARE

- + **User-friendly interface**
Easy to use app and straight-forward results on screen.
- + **Instant reporting**
Scanning results are shown on screen within seconds.
- + **Wireless usage**
NIRLab app pairs with NIRLight via bluetooth and communicates with servers via Wi-Fi or 3G.
- + **Easy and fast download**
The iOS and Android app can be downloaded from Apple or Google store.
- + **Secured cloud**
Complete set of applications communicating with a secured cloud to manage measures and results.



NIRLab Web App

THE SOFTWARE

Desktop app- and browser-accessible platform for data management.

- + **Report history**
Track, manage and compare scans across devices at one place.
- + **Simple data management**
Name, mark, delete or edit analysis results.
- + **Export of data**
Simple data export to Excel sheet.
- + **Save as PDF**
Download the analysis report and save it in PDF format.
- + **User Management**
Organizations and user management tool



PRIVACY

Data Security

- + **Secured data center**
Cloud developed by top-level IT group from the School of Computer Sciences in EPFL, Lausanne, hosted on the university campus secured data center.
- + **Encrypted**
Encrypted communication between mobile app and server.
- + **Full control**
Full control of the information shared in the cloud.
- + **Geolocation**
Geolocation of measurements can be turned on or off.
- + **No sensitive data**
Sample names are coded, and no suspect information is shared.



THE SOLUTION


On-Demand Developments

NIRLAB provides tailor-made services and on-demand development projects to companies across various industries.

- + **Development of new application areas for the identification and quantification of new material and substances**
- + **Integration of new NIR and Raman devices into the existing setting**
- + **White-labelling of our mobile and web application**
- + **Consulting services on data acquisition, data modelling and calibration**


Example:

Falsification of medicine – the case of Viagra



Example:

Fat content in fish



Root mean squared error reduced by 50% from 2.23 to 1.14


How it works



THE APP

How it works

- 1. Connect**
Pairing of NIRLAB app and device is done automatically via Bluetooth connection following two steps:
 - a) Turn on the device
 - b) Open NIRLab app on mobile phone
- 2. Scan**
To perform a scan, point the device on a questioned substance and press the multifunctional button. Scan can be performed with direct contact or through a thin plastic bag.
- 3. Read**
After a few seconds, result of the scan is shown on the screen of your mobile phone.




THE APP SETUP

Calibration

At every start of the app, calibration needs to be performed.

1. To perform a calibration, apply the white reference mirror to the device.
2. Then click on *Calibrate* in the main menu of the app and push the multipurpose button on the device.
3. The process takes a few seconds and is done automatically.

TIP: We recommend to calibrate the device regularly according to the app's notification.




SCANNING MODES

Quick Scan

Click on *Quick Scan* for rapid identification and quantification of a substance.

Procedure:

1. Push the device button
2. Wait a few seconds
3. See the result on screen!



SCANNING MODES

Sample Scan

Click on *Sample Scan* to get an average result of multiples scan to improve accuracy of quantification.

Procedure:

1. Click on *Sample Scan*
2. Fill in information about your sample
3. Push the device button
4. Wait a few seconds

Repeat step 3. and 4. as many times as you need

5. See the averaged result on the screen.



RECOMMENDATION


Scanning Tips

- + **All substances can be measured in direct contact or through thin plastic.** Direct contact produces the best accuracy, especially for quantification.
- + **The sapphire glass should be cleaned** before each scan. This can easily be done by a bit of ethanol on a tissue.
- + **To better assess homogeneity,** *Sample Scans* mode is recommended for powders and high quantity samples.
- + **Small samples** should be measured in an aluminum cup which has a neutral effect on the spectrum.
- + **Point device downwards** when scanning. It is NOT recommended to measure with the device pointing upwards.

TUTORIAL

www.nirlab.com

- [CLICK HERE to watch a tutorial video.](#)



Thank you for your attention

NIRLAB AG / NIRLAB Forensics GmbH

Orsières, Switzerland

+41 21 692 46 57

contact@nirlab.com

www.nirlab.com