



## LED Homecare That Touches the Skin and the Soul A Bold Leap Forward with Yunique H

There's an undeniable thrill in breaking boundaries—and *Yunique H* is leading the charge. This next-generation LED home beauty device is more than just technology; it's a statement.

<https://www.dosinternational.co.kr/>



### The Dawn of Personal Beauty Care in Our Relentless Daily Grind

We're busy. *Too* busy. The so-called “God-saeng” lifestyle—a buzzword among Korea's MZ generation (Millennials + Gen Z)—has become a full-blown cultural trend. People now live by the minute, even by the second, optimizing every hour of the day. The “miracle morning” routine is the new normal, and social media is flooded with productivity challenges and time management hacks. Amid this hyper-efficient chaos, who has time for skincare? Modern life feels merciless at times.

Enter: beauty devices—the ultimate secret weapon. These innovative tools have emerged as a lifesaver for the 20s and 30s crowd, offering high-performance skincare at home without the time commitment or salon price tag. No appointments, no commutes—just results.

In fact, the domestic beauty device market in Korea has more than tripled in size between 2018 and 2022, riding a wave of explosive growth. Industry experts project the market to continue expanding at an annual rate of over 10% through 2030.

Amid this momentum, a new contender has stepped into the spotlight: *Yunique H* by DOS & YOUNIQUE. If you've ever experienced a professional facial, the sleek LED shape of the Yunique H may look familiar. But while traditional LED machines are bulky and expensive—hardly suitable for home use—Yunique H has been engineered to address that very gap.

With a compact form factor and accessible pricing, the Yunique H marks a fresh chapter in at-home LED skincare—making professional-level care possible without stepping outside your door.

## Effortless Skincare at Home LED Therapy Meets Homecare with Younique H

When we think of at-home LED skincare, the classic face mask shape usually comes to mind. The concept is simple yet powerful: targeted wavelengths of light are used to stimulate skin at a cellular level—boosting collagen and elastin production, soothing inflammation, and improving overall tone and texture.

*Younique H* follows this same trusted principle but takes it a step further. With the widest spectrum of LED light colors available in Korea, it offers a fully customizable experience tailored to your *Younique* skin type and condition.

What sets it apart? Its touch-free technology gently activates your skin’s natural potential—making it safe even for sensitive skin. And here’s the game-changer: *Younique H* features an ultra-fine nano mist spray that hydrates and revitalizes the skin while you receive your LED treatment.

It’s not just LED therapy—it’s hydration, radiance, and renewal, all in one sleek device.

### POINT 1: Korea’s First with 7 LED Wavelengths Personalized Skincare for Every Concern

*Younique H* stands out with its industry-leading spectrum of seven distinct LED wavelengths—the most offered in any Korean homecare device. From the gentle violet light (360–430nm) to the deeply penetrating red light (up to 760nm), each wavelength targets different layers of the skin to address a wide range of skin concerns.

Users can select the ideal light color based on their *Younique* needs. Among the most powerful is red light (647–760nm), which reaches deep into the dermis to stimulate fibroblasts—the cells responsible for producing collagen, elastin, and hyaluronic acid. The result? Enhanced skin regeneration, improved elasticity, and a more even, radiant tone.

Whether you’re targeting fine lines, dullness, or overall skin resilience, this customizable light therapy system offers a safe and effective solution—regardless of skin type.



Wavelength Color	Wavelength Range	Main Effects	Recommended Skin Type
Red	647~760nm	Enhances elasticity, reduces wrinkles, promotes skin regeneration	Suitable for all skin types
Blue	430~455nm	Improves acne-prone skin, regulates sebum	Suitable for acne-prone and oily skin
Green	492~550nm	Brightening, pigmentation care	Suitable for pigmented or flushed skin
Yellow	550~580nm	Soothing, reduces puffiness	Suitable for sensitive and flushed skin
Orange	588~647nm	Comprehensive care	Suitable for all skin types
Purple	360~430nm	Improves combination skin issues	Suitable for combination skin
Sky Blue	480~490nm	Calming, regulates skin hydration	Suitable for sensitive skin

\* The above co



## **POINT 2: Ultrasonic Nano Mist Technology for Deep, Dewy Hydration** *Micro-Fine Moisture, Instantly Delivered*

One of *Younique H*'s most distinctive features is its built-in nano mist technology, powered by ultrasonic vibration. This function generates ultra-fine mist particles at the nano level, allowing instant hydration that goes far beyond the surface. Unlike conventional sprays, nano-sized droplets penetrate deep into the skin, delivering moisture right where it's needed—leaving your complexion soft, plump, and visibly refreshed.

The mist feature is fully adjustable with three intensity levels and can be used independently or in tandem with the LED light therapy—creating a perfectly tailored skincare experience.

For an extra glow, use the mist after gentle exfoliation to boost absorption. *Younique H* recommends using purified or mineral water for optimal results—avoid highly concentrated toners, essences, or oil-based liquids.

This is hydration—elevated.

## **POINT 3: Amplify the Effects with LED Therapy** **From Face to Hands—Glow, Firmness, and Radiance All in One**

Want to take your skincare results to the next level? Pairing *Younique H* with the right products can make all the difference. Start with the *D&You Collagenics Face Serum*—a high-concentration formula designed specifically for the face. Apply it before your LED session for visibly brighter, more luminous skin.

Looking to enhance your glow even further? Add the *D&You Whitening Sheet Mask* to your routine. After a gentle exfoliation, apply the mask and follow with a 30-minute LED session. The synergistic effect of the LED light and the active ingredients in the mask helps the formula absorb more deeply—leaving your skin smoother, clearer, and beautifully radiant.

And don't stop at your face—*Younique H* works wonders on your hands, too. If your hands are feeling dry or rough, treat them to an LED session. You'll notice softer, more supple skin in no time.

With the right combinations, beauty really can be head to hand.

### **UNIQUE H — Safe & Easy Usage Tips**

Use the device for about 30 minutes per session. You can adjust the time in 5-minute steps, up to 60 minutes max.

For best results, use it daily. But if you're busy, 2–3 times a week is enough to see benefits. For better results, try using it with sheet masks, serums, or lifting kits to boost absorption and glow.

**Important:** Always wear protective goggles when using the LED light to keep your eyes safe.